# ARCEOLE

CE 0459

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Set for infusion of ophthalmic gas

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**ARCEOLE • INDICATIONS :** Endotamponade of the retina.

### DESCRIPTION:

ARCEOLE is composed of a sterile unidose cylinder of 30ml non sterile gas for ophthalmic surgery and a sterile plastic graduated 50ml syringe mounted with a sterilizing 0.22µm filter and a device to connect it onto the cylinder. A 30G needle is supplied for the intraocular injection. A patient wrist band allows the identification of intraocular gas wearer. A set of 5 patient labels is provided to allow traceability.

The different accessories provided with ARCEOLE have been validated for the herunder detailed use, it is not advisable to substitute them by others.

### COMPOSITION :

ARCEOLE is a sterile unidose cylinder of one of the following gases : SF<sub>6</sub> (Sulfur hexafluoride), C<sub>2</sub>F<sub>6</sub> (Hexafluorohexane), C<sub>3</sub>F<sub>8</sub> (Octafluoropropane).

- Each sterile unidose cylinder contains a pure non sterile gas in a sterile packaging.
- Nb : The gases included in ARCEOLE are non toxic, inert, non flammable, colorless and odorless compounds. Once implanted, they are not metabolized and gradually eliminated through bloodstream and the respiratory tract.
- Expansion and reabsorption rates based on bibliography data :

Gas	Expansion rate	Time allowed for expansion in days	Average tamponade duration in days	Non expansive concentration in %
SF <sub>6</sub>	1.9-2	2	10-15	20
C <sub>2</sub> F <sub>6</sub>	3.3	3	30-35	16-17
C <sub>3</sub> F <sub>8</sub>	4	3	55-65	12

### PRECAUTIONS AND WARNINGS:

• Although non toxic those gases are competing with oxygen when inhaled. Avoid intentional inhalation.

• **It is imperative to stop the anaesthesia with nitrogen protoxide inhalation at least 15 mn before starting the gas injection into the patient eye.**

• The gases for vitreo retinal surgery are exclusively intended to be used by ophthalmologists trained to this technique.

• ARCEOLE is a medical device designed to inject an ophthalmic gas mixed with air,

In case of injection of pure gas under the physician responsibility,

The bibliography<sup>(1)</sup> recommends not to exceed the following quantity:

SF <sub>6</sub>	C <sub>2</sub> F <sub>6</sub>	C <sub>3</sub> F <sub>8</sub>
1.5 ml	1 ml	0.8 ml

• Avoid injection of pure gas in case of narrow irido-corneal angle, a pre-operative gonioscopy is recommended

<sup>(1)</sup>JACOBS PM, TWOMEY JM, LEAVER PK. Behaviour of intraocular gases. Eye 1988, 2, 660-3

• After implantation, a treatment against ocular hypertension should be established and IOP daily monitored.

• Before and after injection of gas, monitoring of the retina central artery is required

• Perform regular follow up of glaucoma patients or patients with IOP, the elderly patients and those with diabetes.

• As a rule, it is necessary to check the ability of the patient to hold his head in the required position during the treatment.

### COMPLICATIONS AND SIDE EFFECTS:

The complications described in the literature are those linked to vitreo-retinal surgery and particularly:

- sub-retinal passage, haemorrhage, haematoma, choroid tear, sub-conjunctival passage.

- Crystalline lens or retina touch with the needle during the injection.

- Risk of separation or capture of the gas bubble between the crystalline lens and the pars plana.

• Potential side effects are: hypertony, retinal tear, cataract, inflammation, endophtalmie, keratopathy and must be considered by the surgeon with regards to the expected benefits.

**A PATIENT IMPLANTED WITH INTRAOCULAR GAS MUST ABSOLUTELY AVOID PRESSURE VARIATION**

**(elevation, travel by plane, diving with or without isobar caisson, hyperbar caisson treatment)**

**A PATIENT IMPLANTED WITH INTRAOCULAR GAS MUST ABSOLUTELY AVOID NITROGEN PROTOXIDE ANAESTHESIA**

In order to remind those side effects and warnings the patient is encouraged to wear the wrist band included in the package.

### INSTRUCTIONS FOR USE:

WARNING the gas is not sterile and must be sterilized by passing through the sterilizing filter supplied with the product. Preparation must be performed extemporaneously.

• Position the syringe equipped with the filter and connector onto the protective seal of the valve and push vigorously the tip of the connector through the seal inside the valve.

• The gas is released inside the syringe and pushes back the plunger, it is sterilized by passing through the 0.22µm filter.

• Disconnect the syringe-filter part from the connector holding the cylinder and swinging the syringe with a right-left movement. (The filter and connector must stay connected to the syringe). Do not pull to disconnect.

• Adjust the plunger to the wished volume of pure gas, using the syringe graduation.

The table below provides help to evaluate the gas concentration:

Gas	Non expansive concentration in %	Volume of pure gas to fill a 50 ml syringe	Volume of air
SF <sub>6</sub>	20	10 ml	40 ml
C <sub>2</sub> F <sub>6</sub>	16	8 ml	42 ml
C <sub>3</sub> F <sub>8</sub>	12	6 ml	44 ml

• Take air in the syringe to obtain the desired proportion mix. The air is sterilized by passing through the syringe filter,

• Remove the filter from the syringe,

• Immediately mount a sterile injection canula with its protective cap, the mix is ready for intraocular injection.

### STORAGE CONDITIONS:

Avoid high temperatures above 50°C/112°F

Store in dry, aerated area.

Do not use if the sterile pouch has been compromised.

Do not use after expiration date.

Single use

Do not re-sterilize

### STERILIZATION:

NOTE: The gas contained in ARCEOLE is not sterile, only the cylinder and accessories are sterilized with ethylene oxide.